# Medtronic

---

## 6726

DF-1 Y-Adaptor/Extender Kit

Kit adaptateur/prolongateur DF-1 Y

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Technical manual • Manuel technique

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Manual técnico

! USA

**Caution:** Federal Law  
(USA) restricts this  
device to sale by or  
on the order of a  
physician.

**CE**  
**0123**  
2001

The following are trademarks or registered trademarks of Medtronic in the United States and possibly in other countries:

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


Medtronic

*Sterile package opening instructions*

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Conformité Européenne  
(European Conformity).

This symbol means that  
the device fully complies  
with European Directive  
AIMD 90/385/EEC.

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

## **Description**

The Medtronic Model 6726 DF-1<sup>1</sup> Y-Adaptor/Extender Kit is designed to adapt two DF-1 lead connectors to a single DF-1 compatible connector port and/or extend the length of a DF-1 lead connector.

Setscrews for securing the lead connector pins are retained within the connector blocks of the adaptor/extender. Grommets seal and insulate the setscrews.

## **Contents of package**

- 1 Model 6726 DF-1 Y-Adaptor/Extender
- 1 torque wrench
- 1 DF-1 pin plug
- Product literature

<sup>1</sup> DF-1 refers to the International Connector Standard (ISO 11318), whereby pulse generators and leads so designated are assured of a basic mechanical fit.

## **Intended use**

The Model 6726 DF-1 Y-Adaptor/Extender is intended to adapt two DF-1 lead connectors to a single DF-1 compatible connector port.

The Model 6726 DF-1 Y-Adaptor/Extender is intended for use when additional extension of a DF-1 lead connector is needed.

The Model 6726 DF-1 Y-Adaptor/Extender kit is intended for single use only.

## **Contraindications**

To date, there are no known contraindications to the use of the Model 6726 DF-1 Y-Adaptor/Extender.

## **Warnings**

An implanted lead or a lead adaptor/extender forms a direct low-resistance current path to the myocardium. During lead implantation and testing, use only battery-powered equipment or line-powered equipment specifically designed for this purpose to protect against fibrillation that may be caused by alternating currents. Line-powered

equipment used in the vicinity of the patient must be properly grounded. Lead connector pins must be insulated from any leakage currents that may arise from nearby line-powered equipment.

Because the lead adaptor/extender is a functional extension of the lead, it is recommended that the lead manual be reviewed for all appropriate warnings, complications, and precautions.

## **Precautions**

### **Non-Medtronic Devices**

Medtronic does not recommend the use of the Model 6726 with any device other than a commercially available implantable pulse generator or implantable cardioverter defibrillator with which it has been tested and demonstrated to be safe and effective.

### **Necessary hospital equipment**

Keep external defibrillation equipment nearby for immediate use during the acute lead system testing, implantation procedure or whenever arrhythmias are possible or intentionally induced during post-implant testing.

## **Handling the lead adaptor/ extender**

Lead adaptor/extender insulators attract small particles such as lint and dust; therefore, protect the lead adaptor/extender from materials shedding these substances.

## **Ensuring product integrity**

Medtronic has sterilized the package contents with ethylene oxide before shipment. This device is for single use only and is not intended to be resterilized.

Carefully inspect the package before opening. If the seal or package is damaged, contact your local Medtronic representative.

Do not use the product after its expiration date.

## **Potential complications**

The potential complications related to the use of a lead system, of which a lead adaptor/extender is a part, include, but are not limited to, the following patient-related conditions: infection, hemothorax, and bleeding after lead placement, after repositioning, or if sutures tear out.

Other potential complications related to the lead system and the programmed parameters of an implantable device include, but are not limited to, the following:

<b>Complications</b>	<b>Possible effects</b>
Displacement, lead system insulation failure or conductor fracture, poor connection to the implantable device or lead.	Intermittent or continuous loss of defibrillation or cardioversion therapy; possible muscle stimulation in the pocket area.
Electrode fracture	Intermittent or continuous loss of defibrillation or cardioversion therapy.
Threshold elevation	Loss of defibrillation or cardioversion therapy.

## **Instructions for use**

Before using the lead adaptor/extender, detach the lead from the implantable device and follow the procedure listed below:

1. **Caution:** Disable the therapy features of the implantable device while the lead is detached from the device. If the therapy features is not disabled while the lead is detached,

there may be an intermittent connection which could result in inappropriate therapy being delivered to the patient. Turn the therapy features back on after connecting the adaptor/extender to the connector block of the device.


Examine the lead for signs of connector pin corrosion, pitting, or any other visible deterioration and test the lead for adequate electrical performance.

**Caution:** If the existing lead is in any way unreliable, do not use the lead. Chronic removal or repositioning of leads may be difficult because of fibrotic tissue development. If a lead must be removed or repositioned, proceed with extreme caution.

- Lead removal may result in avulsion of the endocardium, valve, or vein.
- Lead junctions may separate, leaving the lead tip and bare wire in the heart or vein.
- Cap abandoned leads to avoid transmitting electrical signals.
- For leads that have been severed, seal the remaining lead end and suture the lead to adjacent tissue.


2. Gently insert the connector end of the high voltage lead into either port of the adaptor/extender until the lead connector pin is fully seated in the port (Figure 1). If necessary, use sterile water as an insertion lubricant.

**Figure 1.**

3. Insert the torque wrench through the grommets (Figure 2). Use the torque wrench packaged with the adaptor/extender to turn the setscrews, which are retained in the adaptor/extender connector blocks, until the lead connector pin is firmly held in place and the torque wrench clicks once.

**Figure 2.**

**Note:** If the adaptor/extender is being used primarily as an extender, insert the DF-1 pin plug packaged with the adaptor/extender into either port of the adaptor/extender before inserting the lead into the other port. Tighten the setscrews on the DF-1 pin plug until the pin plug is firmly held in place and the torque wrench clicks once.

4. Gently pull on the lead and adaptor/extender to test the strength of the connection.
5. Test the integrity of the connection by taking defibrillation efficacy measurements according to the directions in the lead product literature.
6. Insert the connector end of the adaptor/extender into the DF-1 port of the implantable device and tighten the setscrews as described in the implantable device product literature.
7. Loosely wrap any excess lead and adaptor/extender length around the implantable device. Insert the lead, adaptor/extender and implantable device into the pocket as indicated in the lead product literature.

**Caution:** Ensure that the adaptor/ extender does not leave the implantable device at an acute angle.

## Specifications (Nominal)

Parameter	Model 6726
Type	Unipolar
Connector	DF-1
Length	15 to 65 cm
Materials	
• Connector Pin	Stainless Steel
• Insulator	Silicone Rubber
• Conductor Cable	MP35N/Ag Alloy
DC Resistance	$\leq 2$ ohms (65 cm)

## Special notice

Medtronic lead adaptor/extender kits consist of lead adaptor/extendors and tools to connect the adaptor/extender to implantable leads and implantable cardioverter defibrillators. Adaptor/ extendors are used with leads, which are implanted in the extremely hostile environment of the human body. Leads, when used with adaptor/extendors, may not be the functional equivalent of leads

that are used without adaptor/extendors. Leads or adaptor/extendors may fail to function for a variety of causes, including, but not limited to: medical complications, body rejection phenomena, allergic reaction, fibrotic tissue, or failure of leads by breakage or by breach of their insulation covering. In addition, adaptor/extendors and tools may easily be damaged by improper handling or use. For leads used with adaptor/extendors, no representation or warranty is made that failure or cessation of function will not occur, or that the body will not react adversely to the implantation of leads or that medical complications (including perforation of the heart) will not occur.

## **Medtronic disclaimer of warranty**

For more information, see the card enclosed in this package.

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© Medtronic, Inc. 2010
Medtronic B.V.
Medtronic Australasia Pty Ltd